

Water Bath

NSW-130 TISSUE FLOTATION BATH

CONSTRUCTION:

- Double walled. Inside made of Aluminium.
- Outside of Mild Steel sheet painted with Epoxy Powder Coating.
- Thermostatic control from few °C above ambient temp. to 80°C ± 1°C.
- To work on 220/230V AC
- Supplied without thermometer.



NSW-130

SIZE INSIDE CHAMBER:

MODEL	SIZE (MM)		LTRS	RATING
	DIA	DEPTH		
TFB-1	210 x 75 mm		3	500 W
TFB-2	250 x 125 mm		6	750 W
TFB-3 New Model	210 x 50 mm		2	500 W
TFB-4 New Model	250 x 50 mm		2	750 W

OPTIONAL:

- Microprocessor based controller cum Indicator with LED display instead of thermostat to cost extra



NSW-130 NEW MODEL