

# Tray Dryer / Drying Oven

## NSW-148

### TRAY DRYER (FOR SEED, LEAVES, VEGETABLES & FRUIT DRYING)

**CONSTRUCTION:**

- This is triple walled, inside made of stainless steel 304 qty and outside made of mild steel outside painted in dark chocolate colour.
- Insulation between the two walls is filled with mineral wool.
- Fitted with Microprocessed electronic digital temperature indicator cum controller with LED display.
- Temp. Range 10°C. Above ambient to 250°C ± 3°C.
- With circulation fan motor and one no. M.S. Drip tray.
- The door has Silicone / Neoprene rubber gasket.



NSW-148

**SIZE INSIDE CHAMBER (STAINLESS STEEL 304 QLTY)**

SIZE (MM)			LTRS	NO OF TRAYS (MADE OF SS 304 QLTY)	RATING
W	D	H			
455	455	455	95	3	2.5 KW
455	455	605	125	4	2.5 KW
605	605	605	224	6	4.0 KW
605	605	910	336	9	4.0 KW

**OPTIONAL:**

- Microprocessed PID Controller having 4 ramp & 4 soak profile having set mode, alarm parameter, configuration parameter, parameter, set point profile, Pv error indicator, parameter for locking etc. to cost extra.
- Alarm: It is possible to provide electronic alarm if, temperature deviates more than ± 5°C to cost extra.
- Automatic timer electronic 999 mts. to cost extra.
- Data logger 4/8 channel temperature recording with basic software and RS-485 PC communication, data can be collect through USB Port to cost extra.