

Tray Dryer / Drying Oven

NSW-149

DRYING OVEN

(FOR SEED, LEAVES, VEGETABLES & FRUIT, MUSHROOM ETC.)

CONSTRUCTION:

- This is triple walled, inside and outside mild steel painted in dark chocolate colour
- Insulation between the two walls is filled with mineral wool.
- It is fitted with Microprocessed Electronic Digital temperature indicator cum controller.
- Circulation fan also to get uniform temperature form plus/minus 3°C. above ambient to 250°C.
- Three Nos. Heavy duty wire mesh mild steel shelves.
- To work on 220/440V AC Supply



NSW-149

SIZE OF DRYER INSIDE CHAMBER:

MODEL	CU. FT	SIZE (MM)			LTRS
		W	D	H	
DOMS-1	18 Cu. Ft	900	600	900	510
DOMS-2	27 Cu. Ft	900	900	900	765
DOMS-3	36 Cu. Ft	900	900	1200	1020
DOMS-4	49 Cu. Ft	900	900	1500	1270
DOMS-5	64 Cu. Ft	1200	1200	1200	1810
DOMS-6	96 Cu. Ft	1200	1200	1800	2720

OPTIONAL:

- Microprocessed PID Controller having 4 ramp & 4 soak profile having set mode, alarm parameter, configuration parameter, profile parameter, set point profile, PV error indicator, parameter for locking etc. to cost extra
- Alarm: It is possible to provide electronic alarm if, temperature deviates more than $\pm 5^{\circ}\text{C}$. to cost extra
- Automatic timer electronic 999 mts. To extra cost
- Electronic Hooter (three tone) to cost extra.
- Automatic timer electronic 9999 mts to cost extra.
- **Extra price for Inner chamber & Outer Chamber SS-304 Qlty to cost extra.**