

# Freezer

## NSW-180 BLOOD BANK REFRIGERATOR "CALTAN"

Outer body and inside chamber & trays are made of stainless steel 304 qty. It has microprocessed digital temp. Controller to control temperature at  $4^{\circ}\text{C} \pm 1^{\circ}\text{C}$ . It will have alarm if door is opened, temperature deviations to low or high. It may be supplied with pull out drawers. It has hermetic compressor with non CFC refrigerant with chart recorder and with automatic voltage stabilizer.

**To works on 220V AC, 50 Hz.**



NSW-180

### VERTICAL MODEL

MODEL	NO. OF BAGS	BLOOD BAGS CAPACITY	NO OF DRAWER/TRAYS
BBR 1 – 4 Cu Ft.	Capacity 60	350ml.	1
BBR 2 – 6 Cu Ft.	Capacity 120	350ml.	2
BBR 3 – 10 Cu Ft.	Capacity 240	350ml.	4
BBR 4 – 12 Cu Ft.	Capacity 360	350ml.	6
BBR 5 – 15 Cu Ft.	Capacity 430	350ml.	6

## NSW-181 PLASMA FREEZER "CALTAN"

Outer made of Mild Steel and inside made of stainless steel 304 qty. Insulation between the two walls are done by PUF insulation. Temperature from ambient to  $-40^{\circ}\text{C}$  is be controlled by microprocessed digital temperature controller. It has low & high temperature cut in the range of  $1^{\circ}\text{C}$  fitted with hermetic compressor with CFC free refrigeration system. It is supplied with automatic voltage stabilizer & chart recorder.

**To work on 220V AC 50 Hz**

### VERTICAL MODEL

MODEL	SIZE (MM)		
	W	D	H
PLFR 1 – 4 Cu Ft.	455	410	610
PLFR 2 – 6 Cu Ft.	505	415	830
PLFR 3 – 10 Cu Ft.	570	550	875
PLFR 4 – 12 Cu Ft.	650	580	900
PLFR 5 – 15 Cu Ft.	700	640	900



NSW-181