

Shakers

NSW- 200 ROTARY SHAKER (VARIABLE SPEED) “ ORBITAL SHAKER”

USES

- These shakers are generally required in the life sciences & Biotechnology
- This series of shaker include Hybridization Incubator.
- The laboratory series of shakers are designed critically of operation for the typical research job

CONSTRUCTION:

- This Shaker's can be supplied in orbital of 3/8", fixed or Variable speed & without RPM Indicator.
- Without flask holding clamps.
- With suitable rating geared motor and electronic speed controller
- These shakers have speed range of 50-350 R.P.M.
- These shakers are complete with adjustable rods from two sides on platform to keep conical flask of any size.
- Which could be adjustable rods according to requirement.
- In built mechanical timer with adjustable 0 to 60 minutes timer with arrangement for continuous operation.
- To work on 220/230V AC.



NSW-200



NSW-200 DIGITAL

SIZE OF SHAKER PLATFORM WITH VARIABLE SPEED

| MODEL | SIZE (MM) W x L | RATING |
|-------|--------------------|--------|
| VS-1 | 380 x 380 | 250 W |
| VS-2 | 455 x 455 | 250 W |
| VS-3 | 605 x 605 | 250 W |
| VS-4 | 910 x 910 | 250 W |
| VS-5 | 1200 x 1200 | 500 W |

OPTIONAL ARRANGEMENTS:

- Digital RPM METER to cost extra
- Digital Timer 999 mts.to cost extra