

Shakers

NSW-200 ROTATOR SHAKER (FIXED SPEED) “VRDL SHAKER”

USES

- Blood bank
- Hematology laboratories
- V.D.R.L. Test and other test.

Operates at 130 or 180 R.P.M. (As desired) Rotates specimen in horizontal plane in 3/4” circle. The 30 x 30 cm platform can accommodate slides, blood bottles, flasks and beakers by use of spring bottle holder furnished with each unit.

CONSTRUCTION:

- Durable all steel construction, Body painted in epoxy powder coating.
- Platform rigidly supported by steel.
- Mounted on four rubber suction feet to prevent creeping.
- 1/35 H.P. Motor
- Without flask holding clamps.
- Operates on 220/230V AC



NSW-200

SIZE OF SHAKER:

MODEL	SIZE (MM) W x L
FSS-1	300 x 300

OPTIONAL:

- The shaker can be supplied with adjustable mechanical 0 to 30/60 minutes or 99 minutes digital timer with arrangement for continuous operation at an extra cost.

NOTE: SHAKING PLATFORM TO ACCOMODATE OTHER CAPACITY OF FLASK

SIZE OF TRAY	FLASK	50 ML	100 ML	250 ML	500 ML	1000 ML	2000 ML
380 x 380 MM	TO HOLD	25 NOS	16 NOS	9 NOS	5 NOS	4 NOS	4 NOS
455 x 455 MM	TO HOLD	38 NOS	25 NOS	16 NOS	9 NOS	5 NOS	5 NOS
605 x 605 MM	TO HOLD	64 NOS	49 NOS	25 NOS	16 NOS	9 NOS	9 NOS
910 x 910 MM	TO HOLD	144 NOS	100 NOS	64 NOS	38 NOS	25 NOS	25 NOS
1200 x 1200 MM	TO HOLD	289 NOS	186 NOS	121 NOS	81 NOS	49 NOS	49 NOS