

Autoclave

NSW-227

VERTICAL AUTOCLAVE "SMART"

- Other specification shall remain as specified in vertical model

MODEL
VAS-3
VAS-4
VAS-5
VAS-6
VAS-7



NSW-227 (SMART)

NSW-229

STERILIZER PRESSURE COOKER TYPE (ELECTRICALLY OPERATED)

- This sterilizer is made of Aluminium for quick transfer of heat which generates steam under pressure.
- Complete made of Aluminium with Aluminium bowl.
- It has double protection for safety.
- Supplied with steel stand
- Electrically operated./ Non- Electric



NSW-229

MODEL ELECTRICALLY OPERATED	SIZE (MM) DIA DEPTH	LTRS	RATING
SPCT-1	230 x 230	10	2.0 KW
SPCT-2	300 x 300	24	2.0 KW

MODEL NON-ELECTRIC	SIZE (MM) DIA DEPTH	LTRS
SPCT	300 x 300	24