

Autoclave

NSW-227 VERTICAL AUTOCLAVE

SUITABLE FOR

- Sterilization of surgical Instruments,
- Dressing material, linen, rubber, plastic material
- In pressure range of 15 to 22 PSI can be controlled

CONSTRUCTION

- With boiler made of Stainless Steel 16 gauge sheet.
- Outer shall is also made of Stainless Steel Sheet
- Lid is made of Stainless Steel Sheet and tightened around by radial locking.
- Neoprene rubber gasket.
- Fitted with water level arrangement to indicate water position inside the boiler
- Pressure gauge, Air/Steam release cock two safety valve are also provided.
- With pizzostat may be set between 10 to 22 PSI.
- Immersion type heating element
- Cut-off automatically when the autoclave is dry.
- Supplied complete with S.S. basket, cord and plug
- Work on 220V AC, 50 Hz supply.



NSW-227 MANUAL

NSW-227 AUTOMATIC

COMPLETE MADE OF STAINLESS STEEL 304 QLTY.

MODEL	SIZE (MM)		LTRS	RATING
	Dia	Depth		
VA-2	300	500 mm	40	3.0 K.W
VA-3	355	550 mm	55	4.0 K.W
VA-4	400	600 mm	80	4.0 K.W
VA-5	450	600 mm	100	4.0 K.W
VA-6	550	600 mm	147	6.0 K.W
VA-7	550	750 mm	187	6.0 K.W

INSIDE & LID MADE OF STAINLESS STEEL 304 QLTY. OUTER MADE OF MILD STEEL POWDER COATED IN CHOCOLATE COLOUR

MODEL	SIZE (MM)		LTRS	RATING
	Dia	Depth		
MA-1	250	450 mm	23	3.0 K.W
MA-2	300	500 mm	40	3.0 K.W
MA-3	355	550 mm	55	4.0 K.W
MA-4	400	600 mm	80	4.0 K.W
MA-5	450	600 mm	100	4.0 K.W
MA-6	550	600 mm	147	6.0 K.W
MA-7	550	750 mm	187	6.0 K.W

OPTIONAL: Fully automatic by PID Controller having facility for:-

- Auto purging: Set Temp.
- Purging Stop: Set Temp.
- Auto Exhaust: Cool
- Auto cycle stop
- Buzzer At cycle over to cost extra
- Automatic timer electronic 999 mts. to cost extra.
- Stainless steel basket (Perforated) to cost extra.
- Micro PLC controller with Data logger, 8 channel recorder with basic software and RS-485 PC communication, pendrive to cost extra.
- Elapse Time
- Set Temperature: Could be set
- Facility for temp. control