

Freezer

NSW-711 UNIQUE LOW FREEZERS (MINUS 86°C) "CALTAN" HORIZONTAL FREEZER (CHEST CABINET)

- Other specification shall remain as specified in Vertical (Upright Model).
- Heavy Duty counter balanced lid has keyed lid lock built in for security.
- Double seal lid gasket and sub lids eliminate frost build up.



NSW-711 HORIZONTAL MODEL

SIZE INSIDE CHAMBER

MODEL	CU. FT	SIZE (MM)			LTRS	RATING
		L	W	D		
HLFU-1S	3.3 Cu.Ft	585	410	410	95	5.0 KW
HLFU-1M	4.0 Cu.Ft	455	400	605	112	5.0 KW
HLFU-2	6.1 Cu.Ft	830	505	415	171	5.0 KW

OPTIONAL:

- Co2 Back Up system with electronically controlled valve to maintain temperature – 40°C to – 50°C to cost extra
- High Pressure Co2 cylinder of 10 Ltrs. Complete with neck ring, valve and siphon tube with liquid Co2 to cost extra.
- Alarm: It is possible to provide electronic alarm if, temp. deviates more than ± 2°C to cost extra
- Data logger 4/8 channel temperature recording with basic software and RS-485 PC communication, data can be collect through USB Port to cost extra.