

Bio-Safety

NSW-205

BIOSAFETY CABINET, CLASS II A-2

- The class II biological safety cabinet is an open-front, ventilated cabinet.
- The Class II biological safety cabinet will provide personnel, environment and product protection.

According to the performance criteria:

- US Standard ANSI/NSF 49:2002
- European Standard EN 13469:2000

Advantage:

- Convenient to move
- Work area surrounded by negative pressure
- Timer function
- LCD display
- Filter life inquiry
- Motored front window Air Pattern and Protected area

Parameter of BSC-II A2 Biological Safety cabinet:

Biological safety cabinet class II A2 type

- HEPA Filter
- HEPA Exhaust
- Electrical Receptals
- Euro Type – 2
- UV and FL Lamp
- Differential pressure Gauge
- Vacuums Non Certificate.



NSW-205

MODEL	SIZE (MM)		
	W	D	H
BSC-1, A2X	900	600	600
BSC-2, A2X	1200	600	600
BSC-3, A2X	1400	600	600
BSC-4, A2X	1800	600	600

Description:

Biological safety cabinet class II A2 type with NSF 49 Certified extra charges will be applicable.

Biological safety cabinet class II A2 type

- 1060 X 810 X 2326
- Electrical supply: 220V/50Hz
- UV and FL lamp, Differential pressure gauge
- Support stand with Urethane wheel
- HEPA Filter (Supplied) and HEPA FILTER (Exhausted)
- Electrical receptacle: Euro type, 2 ea
- Vacuum valve
- NSF 49 in U.S

Biological safety cabinet class II A2 type

- 1380 X 810 X 2326
- Electrical supply: 220V/50Hz
- UV and FL lamp, Differential pressure gauge
- Support stand with Urethane wheel
- HEPA Filter (Supplied) and HEPA FILTER (Exhausted)
- Electrical receptacle: Euro type, 2 ea
- Vacuum valve
- NSF 49 in U.S

Biological safety cabinet class II A2 type

- 11690 X 810 X 2326
- Electrical supply: 220V/50Hz
- UV and FL lamp, Differential pressure gauge
- Support stand with Urethane wheel
- HEPA Filter (Supplied) and HEPA FILTER (Exhausted)
- Electrical receptacle: Euro type, 2 ea
- Vacuum valve
- In accordance with NSF 49

Biological safety cabinet class II A2 type

- 11995 X 810 X 2326
- Electrical supply: 220V/50Hz
- UV and FL lamp, Differential pressure gauge
- Support stand with Urethane wheel
- HEPA Filter (Supplied) and HEPA FILTER (Exhausted)
- Electrical receptacle: Euro type, 2 ea
- Vacuum valve
- NSF 49 in U.S