

NSW-152-R DOUBLE DOOR REFRIGERATOR

Refrigerator offer high capacity and excellent temperature control at an affordable price. These refrigerators are suitable for Vaccines, Reagents, Cultures, Medicines, Storing blood plasma and blood samples. Temperature controller includes digital temperature display of set value and process value, Forced air circulation for maximum temperature uniformity. Our refrigerator are having compressor (CFC FREE) and CFC free refrigerant with all its protective device

USED IN

Microbiology, Virology, Pathology, Food Research Laboratory, Biotechnology etc.

CONSTRUCTION

- This is double walled White painted interior and pre-coated steel exterior.
- Supplied with 5 nos. of height-adjustable shelves made of Tempered glass.
- With non-CFC, high density, Polyurethane Foam insulation of 75mm thickness
- Mercury-free LED interior lighting with externally mounted internal light switch
- Compressor (CFC FREE) with all its protective devices
- With condensing coil.
- 100 mm ground clearance
- With 3 step Automatic Voltage Stabilizer.
- With electronic digital voltmeter.
- Cabinet is mounted on four heavy duty Castor wheels
- Work on 230V AC 50 Hz



NARANG SCIENTIFIC WORKS PVT. LTD.

E-mail: gemsalesnsw@gmail.com

Phone: 9818868239

Web: www.nsw-india.com

SPECIFICATION

Make	NSW INDIA
MODEL	LPR-5
Capacity	500 Liters
Temperature Range	1-8°C
Temperature resolution	0.1°C
Accuracy/Uniformity	1.6°C
Temperature Recovery	Max. 5 mins after opening of door for 1 min
No of Door	Double door upright
Type	2 No. of Hinged Toughen glass door
Type of Insulation	High Density PUF
Type of Cooling	Direct
Defrost System	Automatic defrost system
Control	• Microprocessor based Electronic Controller with 7 Segment LED display • Control panel shall have heating & temperature (set & actual) indicators, Power main indicator
Display	• 7 Segment LED display • The controller shall have a facility to display the temperature in ° C
Door	Soft touch grip for full length door handles
Interior	White painted
Exterior	Pre-coated steel
Shelves	5 nos. of height-adjustable shelves
Shelf load	40Kg
High Quality Door Seals	Yes
Lockable Outer and Inner Lids	Yes
Refrigerants	Environmentally friendly SNAP & F-GAS Compliant HC refrigerants
Safety features	• Both audible and visual alarms High/Low Temperature, Door Ajar alarm with Icon, service alarm and icon, low

**NARANG SCIENTIFIC WORKS PVT. LTD.****E-mail:** gemsalesnsw@gmail.com**Phone:** 9818868239**Web:** www.nsw-india.com

	battery alarm • incorporate Rodent proof power cord
Noise level	55 dBA or less
Voltage stabilizer	5 KVA voltage stabilizer
1" access port	Yes
Pressure equalization port	Yes
Castor	4 Nos. of Heavy duty Castor wheels
Power supply	230V AC 50 Hz



NARANG SCIENTIFIC WORKS PVT. LTD.

E-mail: gemsalesnsw@gmail.com

Phone: 9818868239

Web: www.nsw-india.com